

Automated Cleveland Open Cup Flash Point Tester

# aco-100



### **Interchangeable Ignition Source**

aco-100 is capable of running both gas and electric ignition. Depending on the situation, ignition source can be changed from gas to electric or vice versa.

#### Skimmer for Surface Skin (Asphalt Version)

For aco-100 asphalt version (aco-100as), a unique skimmer has been added to automatically remove the surface skin that forms on the sample when bitumen or some other samples are tested.

#### 7-inch color LCD touch screen

Equipped with a touch screen for simple, intuitive operation.

#### **Easy operation**

Simply select a pre-registered test mode. No need to set detailed parameters.

# Simple Reliable Robust

## **Specifications**

Test Standards	ASTM D92, ISO 2592, IP 36 JIS K 2265-4, AASHTO T48
Measuring Range	Room temperature to 400 °C
Temperature Sensor	PT-100 in stainless steel sheath
Flash Detector	Double ionization rings
Ignition Source	Gas ignition with automatic lighting or electric ignition, interchangeable
Display	7 inch touchscreen color LCD display
Test Modes	Default modes according to test methods. Custom modes also available.
I/O ports	USB Port x 2 , Ethernet x 1
Size	230 mm x 470 mm x 405 mm (W x D x H)
Weight	18 kg
Power Supply	AC100~120V 1400VA / AC220~240V 1300VA 50/60Hz

We reserve the right of changes without prior notice.



